

The Cyprus Orthodontic Society
and **3M** Unitek invite you to

3M Unitek Self Ligating Workshop

presented by Dr. Lars Christensen

July 3rd 2009

Four Seasons Hotel • Limassol, Cyprus



CLARITY | SL
SELF-LIGATING APPLIANCE SYSTEM



SMARTCLIP™
SELF-LIGATING APPLIANCE SYSTEM

3M Unitek

3M Unitek Self Ligating Workshops

presented by Dr. Lars Christensen

The Self Ligating workshops consist of a one-day interactive course with lectures and typodont exercises.

Dr. Lars Christensen will share his clinical experiences with the group using 3M Unitek Self Ligating Brackets and modern light cure materials.

The information supplied and the hands-on nature of the course will give participants a usable 'tool-box' of tips and ideas to take back to their own clinical environment.

AGENDA

- 09:30 Registration and coffee
- 10:00 Introduction to the 3M Unitek Self Ligating Appliance System including set-ups
Typodont Exercise I - Bracket placement on crowded case
- 11:30 Coffee break
- 11:45 Levelling and alignment with emphasis on wire selection, engagement and disengagement
Typodont Exercise II - Initial archwire insertion
- 13:00 Lunch
- 14:00 Treatment mechanics
Typodont Exercise III - Rectangular wires
- 15:30 Coffee break
- 15:45 Finishing and detailing
Typodont Exercise IV - Finishing, tips and tricks
- 17:00 Close

Participation fee: 50€

If you require any further information on this event please contact Mr. Constantinos Parisis on 0030 6946060826.



Positioning

The Advantage of Accurate Bracket Positioning.

Bracket positioning is an important factor for efficient orthodontic treatment. Many doctors agree that precise positioning may be even more important in a self-ligating world that tends to have longer treatment intervals and faster tooth movement. This is why the Clarity™ SL Self-Ligating Appliance has unique features designed to help you place brackets precisely.



Visual References for Bracket Positioning

- True twin design
- Open slot, self-ligating bracket with metal slot liner
- Mesial and distal Nitinol self-ligating clips can serve as long axis indicators
- Available with color changing APC™ PLUS adhesive



Tactile References for Bracket Positioning

- Tooth-specific base design
- Compound contour base



CLARITY™ | SL

SELF-LIGATING APPLIANCE SYSTEM

Patient Consultation

Offering Patients an Attractive Option.

Treatment starts with patient acceptance. That's why the Clarity™ SL Self-Ligating Appliance is designed to help you tap important patient interests such as aesthetic treatment and convenience without giving up what you want: *an appliance that delivers results.*



Aesthetic Bracket

- Translucent
- No ligatures to stain



Convenience of Self-Ligation

- Opportunity to extend treatment intervals
- Potential for faster treatment
- Fast appointments with less chair time



Hygiene and Comfort

- Smooth, molded, radiused edges
- Self-Ligating, no ligatures to attract plaque
- Clean, open slot, self-ligating bracket design



Consultation Aids

- Clarity SL Bracket Typodont
- Clarity SL Bracket 15x Study Model
- Clarity SL Practice Enhancement Kit

